EN **FR** **DE** **IT** **ES** **TR**

Metafil Bulk Fill

EN

INSTRUCTIONS

IMPORTANT :
READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE.
KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.

Resin-based dental restorative material

Metafil Bulk Fill

FOR DENTIST'S USE ONLY

1. What is Metafil Bulk Fill?

Metafil Bulk Fill is a resin-based light-cure, radiopaque, flowable restorative material. It can be used for a posterior restoration including occlusal surface. The content of inorganic fillers is 47.6 vol % (70.0 wt %). The particle size of the inorganic fillers is between 0.01-8.00 µm.

Indications:

- Direct/Indirect posterior restoration including occlusal surfaces.

Contraindications:

Patients who have a history of severe allergic reactions to this product, methacrylate monomers, or acrylate monomers.

Side effects:

In rare cases, components of Metafil Bulk Fill may lead to sensitization. In such cases, Metafil Bulk Fill should no longer be used.

2. Components

- ①Metafil Bulk Fill----- 2.5mL ----- Class IIa CE⁰¹²³
Shade: Low-Flow Universal, Low-Flow A2
- ②18G Needle Tip ----- 5 (with 1 Needle Cap)----- Class I CE⁰¹²³

3. Composition

Metafil Bulk Fill is a paste consisted of the following compositions. Methacrylate monomers, Acrylate monomers, Inorganic fillers, Aromatic amine, and the others.

4. Precautions

This product contains substances that may cause allergic reactions. Read the following points below.

1) Sensitivity

①Metafil Bulk Fill should not be used by clinicians or on patients who are methacrylate monomer-sensitive. If signs of irritation such as redness appear, stop using it immediately and consult a physician.

2) Precautions

①Avoid contact of Metafil Bulk Fill with soft tissues, skin or eyes. Contaminated skin or mucosa with Metafil Bulk Fill should be wiped off immediately with alcohol and then rinsed thoroughly with running water.

②If Metafil Bulk Fill enters an eye, immediately rinse the eye thoroughly with running water. The patients should then be examined by an ophthalmologist.

③Avoid accidental ingestion while applying this product. If ingested, the patient must seek medical attention.

④Dentist should wear rubber or PVC dental gloves and eyewear.

⑤Proper use of eye protection is recommended when using light irradiation devices.

⑥Use a dental rubber dam to isolate the treatment area from contamination.

⑦In the case of pulp exposure or deep cavities with close proximity to pulp, cap the pulp with calcium hydroxide etc. Do not use eugenol-containing material.

⑧Needle Tip is for one-time use only. It is not autoclavable or chemically sterilizable. Dispose it after each patient to avoid cross-contamination. After use, discard any excess.

⑨Use the Needle Tips only for this product.

⑩Do not mix the components with other materials. Once dispensed, do not return this product to the containers. Use Metafil Bulk Fill only for the applications recommended in this leaflet.

⑪Product life time can be shortened depending on the circumstances of use. Read instructions of relevant products carefully before use.

5. Storage

- ①Store Metafil Bulk Fill at a room temperature between 1 – 30°C. Avoid high temperatures above 30 °C, high humidity and direct sunlight.
- ②Shelf life: 3 years in a room temperature.
- ③Do not use this product after the expiration date. See outer package for expiration date.

6. Disposal

Dispose of empty containers, packaging, etc. with the normal household waste.

As in any dental treatment, the patient's individual constitution and the unique requirements of each clinical case must be considered before selecting materials and conditions for use.

7. How to use Metafil Bulk Fill

<Metafil Bulk Fill syringe>

①Remove the black cap by twisting it counter-clockwise. Save the black cap.

②Install the Needle Tip into the syringe certainly and make sure the Needle Tip and the syringe are connected tightly. Then, cover the tip with the Needle Cap until dispensing Metafil Bulk Fill.

<Direct restoration>

①Prepare the cavity by removing caries and any existing restorations as indicated.

②Protect pulp by covering the dentin close to the pulp with adequate pulp capping materials.

③Adhesive Bonding System :

Apply light-cured bonding system according to the manufacturer’s instructions. We recommend using Hybrid Bond made by SUN MEDICAL CO., LTD.

④Dispense Metafil Bulk Fill directly into the cavity by slowly squeezing the plunger. Cavities can be filled up to 4mm in one increment.

⑤Light cure Metafil Bulk Fill. Refer to Table 1 for the curing time.

Note: For deeper preparations, add more Metafil Bulk Fill by increments.

⑥Check occlusion with a thin articulating paper. Examine centric and lateral excursion contacts. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.

Table 1: Curing time and curing depth of Metafil Bulk Fill				
Type	Light intensity (mW/cm²)	Curing depth (mm)	Curing Time (sec)	
			Low-Flow Universal	Low-Flow A2
LED	1,000 <	4	10	10
Fast halogen	1,000 <		10	20
Conventional halogen	400 <		20	40
Plasma arc	1,840 <		6	9

⑦Metafil Bulk Fill is polished quickly and easily to a high luster. Most resin discs or polishing disc, paste and polishing systems are acceptable.

When Metafil Bulk Fill is used as a base/liner, please follow the directions below.

<Base/liner under direct restorations>

①Prepare the tooth in the same way as for direct restoration, step ①–⑤.

②Place a surface layer of a methacrylate-based universal/posterior composite resin restorative. When Metafil Bulk Fill is used as a base/liner, surface layer composite resin material should have at least 2 mm thickness.

③For working and finishing directions, follow the manufacturer’s instruction for use.

<Indirect restoration with prosthesis>

①Prepare the cavity by following the step ①–⑤ for direct restoration to place Metafil Bulk Fill as base/liner. Finish the preparation.

②Take an impression with agar/alginate or silicone materials.

③While waiting for the fabrication of final prosthesis, properly protect the tooth surface with temporary restorative material or temporary crown. When using resin-based temporary restorative material, apply dental separating agent to the cavity wall.

④Bond the prepared prosthesis with commercially available dental cement such as Super-Bond C&B or SEcure made by SUN MEDICAL CO., LTD. Prime the prosthesis if indicated by the instruction of the preferred cement. We recommend V-PRIMER for priming precious metal and Super-Bond Universal Ceramic Primer for ceramic/zirconia prosthesis.

Note: Do not use eugenol-containing temporary restorative material.

8. To get the best results with Metafil Bulk Fill

①Use Metafil Bulk Fill with 18G Needle Tip.

②Allow this product to reach room temperature prior to use for proper flow if Metafil Bulk Fill is stored in a cool place.

③Metafil Bulk Fill is light sensitive material. Protect from ambient or operatory lights. Reduce/turn off the light while the restorative procedure as indicated to maintain proper working time.

④If the cavity is very deep, incremental restoration may be required (refer to Table 1.)

⑤Check the performance of the light curing device. Clean and exchange the light bulb periodically. Insufficient polymerization from the compromised irradiator may cause the lower adhesive strength.

⑥All containers and the syringe of Metafil Bulk Fill should be closed tightly immediately after each dosing.

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