




Universal  
scaler tip




Visual wear  
indicator


The Scorpion insert tip is coated with a 2 µm thick and strong titanium nitride layer. To maintain efficiency and for patient comfort the insert has to be changed when the titanium coating has disappeared

Very thin


Very powerful




Scorpion tip



Other tip



Scorpion tip



Other tip

Scorpion tip reaches interproximal spaces without difficulty

A strong cavitation is guaranteed even with a low power scaling device

CLiP  
& INSERT-i


Safe Implant Cleaner

Clean implants and restorations without changing your tip !



CLiP x 5

+




Insert i

+



Dynamometric key

+



CLIPPER

EMS - Woodpecker



Acteon Satelec - NSK varios - DTE



Sirona Sirosonic



Kavo PiezoLux - SonoSoft - Mectron



Kavo SONICflex  
2000-2003



W&H Synea-Proxeo  
NSK Ti-Max S970 - AS2000



Kavo PiezoLED  
(before 2022)  
Kavo PIEZOsoft



LM - Amdent



# CLiP & INSERT-i

## Safe Implant Cleaner

- PEEK tip, high quality medical polymer

- Allows smooth cleaning of implants and restorations
- Does not damage sensitive surfaces
- Sterilizable and reusable



### Easily attached and removed thanks to the CLIPPER

- Essential accessory for fixing and removing the CLiP securely
- CLiPs wear easily visualised
- Sterilizable - May contain 5 CLiPs

### Exclusively for use on Scorpion-i tips



#### New dual function tip

Can be used without the CLiP for supragingival scaling of teeth or with the CLiP for cleaning implants and restorations

Visual wear indicator - Efficiency - Long life

Especially designed and patented to secure the CLiP

Available as a complete kit\*



CLiP x 5

+



i-tip

+



Dynamometric key

+



CLIPPER

Except **I-AIR** and **I-LM** kits: contain 10 CLiPs + i-tip + CLIPPER

\* Or separately (CLiPs available in packs of 5)



4 rue Enrico Fermi - 63540 Romagnat - FRANCE  
E-mail: [contact@scorpion.fr](mailto:contact@scorpion.fr)  
Web: [www.scorpion.fr](http://www.scorpion.fr)

